

**CERTIFICATE OF ANALYSIS**

**Protease-Inhibitor-Mix P**

**Cat.No. : 39103**

**Lot No.: P140399**

<b>Cat. number</b>	<b>Components</b>	<b>Target Protease</b>
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
32010	1,10-Phenanthroline	Metalloproteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
<b>Storage (°C)</b>		-15 to -25
<b>Minimum shelf life</b>		09/ 2017

**Reconstitution** Dissolve the content of the vial (64.2 mg) in 1ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 04.09.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.