

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix B

Cat.No. : 39105

Lot No.: P200279

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
11280	EDTA-Disodium	Metalloproteases
Minimum shelf life		07.08.2024
Storage (°C)		-15 to -25
Reconstitution	Dissolve the content of the vial (26.4 mg) in 1 ml 90 % DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). Slight turbidity in the reconstituted solution is normal. After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).	

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 06.07.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.