

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix HP

Cat.No. : 39106

Lot No.: P120349

| Cat. number | Components | Target Protease |
|--------------------|----------------------------|---|
| 12745 | AEBSF-HCl | Serine Proteases |
| 13718 | Aprotinin from bovine lung | Serine Proteases |
| 21100 | E-64 | Cysteine Proteases |
| 51867 | Leupeptin | Cysteine Proteases and Trypsin-like Serine Proteases |

Storage (°C) -15 to -25

Minimum shelf life 09/ 2015

Reconstitution Dissolve the content of the vial (12.7 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks).

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 04.09.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.