

# CERTIFICATE OF ANALYSIS

**Protease-Inhibitor-Mix HP**

**Cat.No. : 39106**

**Lot No.: P130026**

<b>Cat. number</b>	<b>Components</b>	<b>Target Protease</b>
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases

**Storage (°C)** -15 to -25

**Minimum shelf life** 01/ 2016

**Reconstitution** Dissolve the content of the vial (12.7 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks).

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 28.01.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

This report has been computer generated and does not contain a signature.