

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix HP

Cat.No. : 39106

Lot No.: P160190

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
Storage (°C)		-15 to -25
Minimum shelf life		04.05.2023
Reconstitution	Dissolve the content of the vial (12.7 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks).	

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 11.05.2021

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.