

# CERTIFICATE OF ANALYSIS

**Protease-Inhibitor-Mix HP**

**Cat.No. : 39106**

**Lot No.: P210185**

| <b>Cat. number</b>        | <b>Components</b>   | <b>Target Protease</b>                                  |
|---------------------------|---|---|
| 12745                     | AEBSF-HCl   | Serine Proteases  |
| 13718                     | Aprotinin from bovine lung  | Serine Proteases  |
| 21100                     | E-64  | Cysteine Proteases                                      |
| 51867                     | Leupeptin   | Cysteine Proteases and<br>Trypsin-like Serine Proteases |
| <b>Storage (°C)</b>       |   | -15 to -25  |
| <b>Minimum shelf life</b> |   | 27.05.2024  |
| <b>Reconstitution</b>     | Dissolve the content of the vial (12.7 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks). |   |

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 31.05.2021

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.