

# CERTIFICATE OF ANALYSIS

**Protease-Inhibitor-Mix HP**

**Cat.No. : 39106**

**Lot No.: P230249**

| Cat. number        | Components                                                                                                                                                                                                                        | Target Protease                                         |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| 12745              | AEBSF-HCl                                                                                                                                                                                                                         | Serine Proteases                                        |
| 13718              | Aprotinin from bovine lung                                                                                                                                                                                                        | Serine Proteases                                        |
| 21100              | E-64                                                                                                                                                                                                                              | Cysteine Proteases                                      |
| 51867              | Leupeptin                                                                                                                                                                                                                         | Cysteine Proteases and<br>Trypsin-like Serine Proteases |
| Minimum shelf life |                                                                                                                                                                                                                                   | 30.06.2026                                              |
| Storage (°C)       |                                                                                                                                                                                                                                   | -15 to -25                                              |
| Reconstitution     | Dissolve the content of the vial (12.7 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks). |                                                         |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 05.07.2023**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.