

CERTIFICATE OF ANALYSIS

Protein Test Mixture 4

Cat.No. : 39208

Lot No.: P190111

Protein Components	Molecular weight
Carbonic anhydrase	29 000
Albumin egg	45 000
Albumin bovine	67 000
Phosphorylase B	97 400
Minimum shelf life	13.03.2021
Storage (°C)	+2 to +8
Reconstitution	Dissolve to a concentration of 1mg/ml in SDS-sample buffer. A volume of 5µl is sufficient for SDS-PAGE (in a standard form of 80 x 60 x 0.75mm).

SDS-sample buffer:

125 mM Tris-HCl pH 6.8; 2% SDS; 15% glycerol;
10 mM DTT; 0.025% Bromphenol blue; 0.025% Orange G

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 14.03.2019

Christian Monsler

Marina Konter

This report has been computer generated and does not contain a signature.