

CERTIFICATE OF ANALYSIS



Protein Test Mixture 4

Cat.No. : 39208

Lot No.: P230413

Protein Components	Molecular weight
---------------------------	-------------------------

Carbonic anhydrase	29 000
---------------------------	--------

Albumin egg	45 000
--------------------	--------

Albumin bovine	67 000
-----------------------	--------

Phosphorylase B	97 400
------------------------	--------

Minimum shelf life	22.11.2025
---------------------------	------------

Storage (°C)	+2 to +8
---------------------	----------

Reconstitution	Dissolve to a concentration of 1mg/ml in SDS-sample buffer. A volume of 5µl is sufficient for SDS-PAGE (in a standard form of 80 x 60 x 0.75mm).
-----------------------	---

SDS-sample buffer:

125 mM Tris-HCl pH 6.8; 2% SDS; 15% glycerol;
10 mM DTT; 0.025% Bromphenol blue; 0.025% Orange G

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 24.11.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.