

Certificate of Analysis**Protein Test Mixture 5 for SDS PAGE****Cat.No. : 39209****Lot No.: P100291**

Protein Components	Molecular weight (kDa)
---------------------------	-----------------------------------

Trypsin inhibitor bovine lung	6 500
--	-------

Cytochrom C	12 500
--------------------	--------

Trypsin inhibitor soybean	21 000
--------------------------------------	--------

Carbonic anhydrase	29 000
---------------------------	--------

Minimum shelf life	02.09.2017
---------------------------	------------

Storage (°C)	+2 to +8
---------------------	----------

Reconstitution	Dissolve to a concentration of 1mg/ml in SDS-sample buffer. A volume of 5µl is sufficient for SDS-PAGE (in a standard form of 80 x 60 x 0.75mm).
-----------------------	--

SDS-sample buffer:

125 mM Tris-HCl pH 6.8; 2% SDS; 15% glycerol;
10 mM DTT; 0.025% Bromphenol blue; 0.025% Orange G**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.****Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 16.06.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.