

Certificate of Analysis**Protein Test Mixture 5 for SDS PAGE****Cat.No. : 39209****Lot No.: P170272****Protein Components****Molecular weight****Trypsin inhibitor (bovine lung)**

6 500

Cytochrom C

12 500

Trypsin inhibitor (soybean)

21 000

Carbonic anhydrase

29 000

Minimum shelf life

07.08.2024

Storage (°C)

+2 to +8

Reconstitution

Dissolve to a concentration of 1 mg/ml in SDS-sample buffer.
A volume of 5 µl is sufficient for SDS-PAGE
(in a standard form of 80 x 60 x 0.75mm).

SDS-sample buffer:

125 mM Tris-HCl pH 6.8; 2% SDS; 15% Glycerol;
10 mM DTT; 0.025% Bromphenol blue; 0.025% Orange G

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control**Printing date: 20.06.2023**

Dr. Judith Koch

Marina Konter

This report has been computer generated and does not contain a signature.