

Certificate of Analysis

Protein Test Mixture 5 for SDS PAGE

Cat.No. : 39209

Lot No.: P230331

Protein Components

Molecular weight

Trypsin inhibitor (bovine lung)

6 500

Cytochrom C

12 300

Trypsin inhibitor (soybean)

21 000

Carbonic anhydrase

29 000

Minimum shelf life

07.09.2028

Storage (°C)

+2 to +8

Reconstitution

Dissolve to a concentration of 1 mg/ml in SDS-sample buffer.
A volume of 5 µl is sufficient for SDS-PAGE
(in a standard form of 80 x 60 x 0.75mm).

SDS-sample buffer:

125 mM Tris-HCl pH 6.8; 2% SDS; 15% Glycerol;
10 mM DTT; 0.025% Bromphenol blue; 0.025% Orange G

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 20.09.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.