

Certificate of Analysis

SERVA Native Marker
Liquid Mix

Cat.No. : 39213
Lot No.: P240148

Protein Components	Molecular weight (kDa)
Ferritin (horse)	720 / 450
Phycoerythrin B	240
Lactate dehydrogenase (porcine heart)	146
Albumin (bovine)	67
Albumin (egg)	45
Trypsin inhibitor (soybean)	21
Storage (°C)	-15 to -25
Minimum shelf life	26.04.2025

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 30.04.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.