

# CERTIFICATE OF ANALYSIS

**SERVA Unstained Protein Standard 6.5 - 97 kDa**

**Cat.No. : 39214**

**Lot No. : P210014**

<b>Protein Components</b>	<b>Molecular weight</b>
<b>Phosphorylase B</b>	97 400
<b>Albumin bovine (BSA)</b>	67 000
<b>Albumin egg</b>	45 000
<b>Carbonic anhydrase</b>	29 000
<b>Trypsin inhibitor (soybean)</b>	21 000
<b>Cytochrome C</b>	12 500
<b>Trypsin inhibitor (bovine lung)</b>	6 500
<b>Minimum shelf life</b>	21.01.2022
<b>Storage (°C)</b>	-15 to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 26.01.2021**

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.