

CERTIFICATE OF ANALYSIS

SERVA Unstained Protein Standard 6.5 - 97 kDa

Cat.No. : 39214

Lot No. : P230115

Protein Components

Molecular weight

Phosphorylase B

97 400

Albumin bovine (BSA)

67 000

Albumin egg

45 000

Carbonic anhydrase

29 000

Trypsin inhibitor (soybean)

21 000

Cytochrome C

12 500

Trypsin inhibitor (bovine lung)

6 500

Minimum shelf life

17.03.2024

Storage (°C)

-15 to -25

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 11.01.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.