

**Certificate of Analysis**

---

**SERVA Unstained SDS PAGE Protein Marker**  
6.5 – 200 kDa, Liquid Mix

---

**Cat.No. : 39215**  
**Lot No.: 160948**

<b>Protein Components</b>	<b>Molecular weight (kDa)</b>
<b>Myosin</b> (bovine muscle)	200
<b>b-Galactosidase</b> (E. coli. recombinant)	116
<b>Albumin</b> (bovine)	67
<b>Albumin</b> (egg)	45
<b>Carbonic anhydrase</b> (bovine erythrocytes)	29
<b>Trypsin inhibitor</b> (soybean)	21
<b>Lysozym</b> (egg)	14.3
<b>Aprotinin</b> (bovine lung)	6.5
<b>Storage (°C)</b>	-15 to -25
<b>Minimum shelf life</b>	24.09.2018

<b>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</b>
--

---

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date: 17.10.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.