

Certificate of Analysis**SERVA Unstained SDS PAGE Protein Marker**
6.5 – 200 kDa, Liquid Mix**Cat.No. : 39215**
Lot No.: 180398

Protein Components	Molecular weight (kDa)
Myosin (bovine muscle)	200
b-Galactosidase (E. coli. recombinant)	116
Albumin (bovine)	67
Albumin (egg)	45
Carbonic anhydrase (bovine erythrocytes)	29
Trypsin inhibitor (soybean)	21
Lysozym (egg)	14.3
Aprotinin (bovine lung)	6.5
Storage (°C)	-15 to -25
Minimum shelf life	24.02.2020

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control**Printing date: 11.04.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.