

Certificate of Analysis

SERVA Prestained SDS PAGE Protein Marker
6.5 – 200 kDa, Liquid Mix

Cat.No. : 39216
Lot No.: 170117

Protein Components	Molecular weight (kDa)
Myosin (bovine muscle)	200
b-Galactosidase (E. coli. recombinant)	116
Albumin (bovine)	67
Albumin (egg)	45
Carbonic anhydrase (bovine erythrocytes)	29
Trypsin inhibitor (soybean)	21
Lysozym (egg)	14.3
Aprotinin (bovine lung)	6.5
Storage (°C)	-15 to -25
Minimum shelf life	25.11.2018

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 10.02.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.