

**Certificate of Analysis**

<b>SERVA Native Marker</b>	<b>Cat.No. :</b>	<b>39219</b>
Liquid Mix	<b>Lot No.:</b>	<b>P150445</b>

<b>Protein Components</b>	<b>Molecular weight (kDa)</b>
<b>Ferritin</b> (horse)	720 / 450
<b>Urease</b> (jack bean)	545 / 272
<b>Lactate dehydrogenase</b> (porcine heart)	146
<b>Albumin</b> (bovine)	67
<b>Albumin</b> (egg)	45
<b>Trypsin inhibitor</b> (soybean)	21
<b>Storage</b> (°C)	-15 to -25
<b>Minimum shelf life</b>	26.02.2017

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 14.06.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.