

Certificate of Analysis

| | | |
|----------------------------|------------------|----------------|
| SERVA Native Marker | Cat.No. : | 39219 |
| Liquid Mix | Lot No.: | P170274 |

| Protein Components | Molecular weight (kDa) |
|--|-----------------------------------|
| Ferritin (horse) | 720 / 450 |
| Urease (jack bean) | 545 / 272 |
| Lactate dehydrogenase (porcine heart) | 146 |
| Albumin (bovine) | 67 |
| Albumin (egg) | 45 |
| Trypsin inhibitor (soybean) | 21 |
| Storage (°C) | -15 to -25 |
| Minimum shelf life | 09.11.2018 |

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 15.08.2017

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.