

Certificate of Analysis

SERVA Native Marker
Liquid Mix

Cat.No. : 39219
Lot No.: P210242

Protein Components	Molecular weight (kDa)
Ferritin (horse)	720 / 450
Urease (jack bean)	545 / 272
Lactate dehydrogenase (porcine heart)	146
Albumin (bovine)	67
Albumin (egg)	45
Trypsin inhibitor (soybean)	21
Storage (°C)	-15 to -25
Minimum shelf life	16.10.2022

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 19.07.2021

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.