

CERTIFICATE OF ANALYSIS

SingleQuant Assay Kit

Cat.No. : 39226

Contr.No. : P150302

Kit Components

Storage

Reference standard (BSA)

aliquots at -20°C

Wash solution

room temperature

Elution solution

room temperature

SingleQuant stock solution

+2 to +8 $^{\circ}\text{C}$

Application

protein determination

Minimum shelf life

26.03.2019

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 17.01.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.