

CERTIFICATE OF ANALYSIS

SingleQuant Assay Kit

Cat.No.: 39226
Contr.No.: P230056

Kit Components

Storage

Reference standard (BSA)

aliquots at -20°C

Wash solution

room temperature

Elution solution

room temperature

SingleQuant stock solution

+2 to +8°C

Application

protein determination

Minimum shelf life

22.05.2025

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 27.02.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.