

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

| | |
|------------------------------------|--------------------------------|
| BCA Protein Assay Macro Kit | Kat.Nr./Cat.No. : 39228 |
| BCA Protein Assay Macro Kit | Lot/Contr.No. : P220369 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--------------------------|--|---------------------------|
| Anwendung application | BCA Protein Assay | entspricht corresponds | entspricht corresponds |
| Komponenten components | | Reagent A: BCA-Lösung Reagent B: Kupfersulfat-Lösung Proteinstandard (Rinderserumalbumin) Reagent A: BCA solution Reagent B: copper sulphate solution Protein standard (bovine serum albumin) | entspricht corresponds |
| Mindesthaltbarkeit minimum shelf life | | | 21.04.2024 |
| Lagerung (°C) storage | | | +2 bis/ to +8 |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.11.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.