

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

BCA Protein Assay Micro Kit	Kat.Nr./Cat.No. :	39229
BCA Protein Assay Micro Kit	Lot/Contr.No. :	P170410

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Anwendung application	BCA Protein Assay	entspricht corresponds	entspricht corresponds
Komponenten components		Reagent A: Puffer-Lösung Reagent B: BCA-Lösung Reagent C: Kupfersulfat-Lösung Proteinstandard (Rinderserumalbumin)	entspricht
		Reagent A: Buffer-Lösung Reagent B: BCA solution Reagent C: copper sulphate solution Protein standard (bovine serum albumin)	corresponds
Mindesthaltbarkeit minimum shelf life			11.12.2019
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.11.2018

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.