

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

| | |
|--|--------------------------------|
| SERVA Purple Protein Quantifizierungs-Assay | Kat.Nr./Cat.No. : 39235 |
| SERVA Purple Protein Quantification Assay | Lot/Contr.No. : P190157 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--------------------------|---------------------------------------|---------------------------|
| Aussehen appearance | | gelbe Flüssigkeit yellow liquid | entspricht corresponds |
| Reinheit (%) purity | HPLC | Status | > 95 |
| Mindesthaltbarkeit minimum shelf life | | | 17.08.2023 |
| Lagerung (°C) storage | | | -15 bis/ to -25 |

SERVA Electrophoresis GmbH

| |
|--|
| <p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p> |
|--|

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.08.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.