

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA Purple Protein Quantifizierungs-Assay	Kat.Nr./Cat.No. :	39235
SERVA Purple Protein Quantification Assay	Lot/Contr.No. :	P220391

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbe Flüssigkeit yellow liquid	entspricht corresponds
Reinheit (%) purity	HPLC	Status	> 95
Mindesthaltbarkeit minimum shelf life			28.02.2026
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.07.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.