

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Lowry Assay Kit</b>	<b>Kat.Nr./Cat.No. :</b>	<b>39236</b>
Lowry Assay Kit	<b>Lot/Contr.No. :</b>	<b>P170157</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	Solution A	klare, bläuliche Lösung clear, bluish solution	entspricht corresponds
	Solution B	klare farblose Lösung clear, colorless solution	entspricht corresponds
	Solution C	klare, gelbe Lösung clear, yellow solution	entspricht corresponds
	Albumin bovine Fraction V	perlweißes Lyophilisat oyster white lyophilisate	entspricht corresponds
<b>Anwendung</b> application	Lowry assay	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			15.09.2018
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 19.04.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.