

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Lowry Assay Kit	Kat.Nr./Cat.No. :	39236
Lowry Assay Kit	Lot/Contr.No. :	P190096

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance	Solution A	klare, bläuliche Lösung clear, bluish solution	entspricht corresponds
	Solution B	klare farblose Lösung clear, colorless solution	entspricht corresponds
	Solution C	klare, gelbe Lösung clear, yellow solution	entspricht corresponds
	Albumin bovine Fraction V	perlweißes Lyophilisat oyster white lyophilisate	entspricht corresponds
Anwendung application	Lowry assay	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			27.08.2020
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.03.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.