

**ANALYSENZERTIFIKAT**  
 CERTIFICATE OF ANALYSIS

<b>Lowry Assay Kit</b>	<b>Kat.Nr./Cat.No. : 39236</b>
Lowry Assay Kit	<b>Lot/Contr.No. : P200046</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	Solution A	klare, bläuliche Lösung clear, bluish solution	entspricht corresponds
	Solution B	klare farblose Lösung clear, colorless solution	entspricht corresponds
	Solution C	klare, gelbe Lösung clear, yellow solution	entspricht corresponds
	Albumin bovine Fraction V	perlweißes Lyophilisat oyster white lyophilisate	entspricht corresponds
<b>Anwendung</b> application	Lowry assay	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			21.01.2022
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 23.04.2021**

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.