

CERTIFICATE OF ANALYSIS

SERVA Unstained Protein Standard IV	Cat.No. : 39250
	Lot No. : 171094

Protein Recombinant	Molecular Weight (kDa)
Protein 1	200
Protein 2	150
Protein 3	100
Protein 4	85
Protein 5	60
Protein 6	50
Protein 7	40
Protein 8	30
Protein 9	25
Protein 10	20
Protein 11	15
Protein 12	10
Storage (°C)	-15 to -25
Minimum Shelf Life	18.03.2019

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.12.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.