

**CERTIFICATE OF ANALYSIS**

---

**SERVA Triple Color Protein Standard I**  
15 – 170 kDa, Liquid Mix

---

**Cat.No. : 39251**  
**Lot No. : 250877**

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	170
<b>Protein 2 (blue)</b>	130
<b>Protein 3 (blue)</b>	100
<b>Protein 4 (red)</b>	70
<b>Protein 5 (green)</b>	60
<b>Protein 6 (blue)</b>	45
<b>Protein 7 (blue)</b>	35
<b>Protein 8 (blue)</b>	25
<b>Protein 9 (green)</b>	15
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	06.02.2027

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

---

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date: 30.09.2025**

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.