

**CERTIFICATE OF ANALYSIS****SERVA Dual Color Protein Standard III**  
7 – 240 kDa, Liquid Mix**Cat.No. : 39252**  
**Lot No. : 170696**

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	240
<b>Protein 2 (blue)</b>	140
<b>Protein 3 (blue)</b>	100
<b>Protein 4 (orange)</b>	70
<b>Protein 5 (blue)</b>	50
<b>Protein 6 (blue)</b>	35
<b>Protein 7 (orange)</b>	25
<b>Protein 8 (blue)</b>	20
<b>Protein 9 (blue)</b>	15
<b>Protein 10 (blue)</b>	7
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	19.07.2019

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date: 24.07.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.