

CERTIFICATE OF ANALYSIS

| | |
|--|-------------------------|
| SERVA Dual Color Protein Standard III | Cat.No. : 39252 |
| 7 – 240 kDa, Liquid Mix | Lot No. : 180011 |

| Protein Components | Molecular Weight (kDa) |
|---------------------------|-------------------------------|
| Protein 1 (blue) | 240 |
| Protein 2 (blue) | 140 |
| Protein 3 (blue) | 100 |
| Protein 4 (orange) | 70 |
| Protein 5 (blue) | 50 |
| Protein 6 (blue) | 35 |
| Protein 7 (orange) | 25 |
| Protein 8 (blue) | 20 |
| Protein 9 (blue) | 15 |
| Protein 10 (blue) | 7 |
| Storage (°C) | -15 to -25 |
| Minimum Shelf Life | 04.01.2020 |

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.01.2018

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.