

CERTIFICATE OF ANALYSIS

SERVACHrom Protein Standard III 7 – 240 kDa, Liquid Mix	Cat.No. : 39255 Lot No. : 170697
---	---

Protein Components	Molecular Weight (kDa)
Protein 1 (orange)	240
Protein 2 (blue)	140
Protein 3 (violet)	100
Protein 4 (orange)	70
Protein 5 (yellow)	50
Protein 6 (green)	35
Protein 7 (orange)	25
Protein 8 (blue)	20
Protein 9 (red)	15
Protein 10 (blue)	7
Storage (°C)	-15 to -25
Minimum Shelf Life	19.07.2019

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 24.07.2017

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.