

CERTIFICATE OF ANALYSIS

SERVA Triple Color Protein Standard II	Cat.No. : 39257
11 – 180 kDa, Liquid Mix	Lot No. : 140561

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	180
Protein 2 (blue)	135
Protein 3 (blue)	100
Protein 4 (orange)	75
Protein 5 (blue)	63
Protein 6 (blue)	48
Protein 7 (blue)	35
Protein 8 (green)	25
Protein 9 (blue)	17
Protein 10 (blue)	11
Storage (°C)	-15 to -25
Minimum Shelf Life	05/ 2016

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.05.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.