

**CERTIFICATE OF ANALYSIS**

<b>SERVA Triple Color Protein Standard II</b>	<b>Cat.No. : 39257</b>
11 – 180 kDa, Liquid Mix	<b>Lot No. : 141303</b>

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	180
<b>Protein 2 (blue)</b>	135
<b>Protein 3 (blue)</b>	100
<b>Protein 4 (orange)</b>	75
<b>Protein 5 (blue)</b>	63
<b>Protein 6 (blue)</b>	48
<b>Protein 7 (blue)</b>	35
<b>Protein 8 (green)</b>	25
<b>Protein 9 (blue)</b>	17
<b>Protein 10 (blue)</b>	11
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	12/ 2016

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 09.01.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.