

CERTIFICATE OF ANALYSIS



SERVA Triple Color Protein Standard II
11 – 180 kDa, Liquid Mix

Cat.No. : 39257
Lot No. : 240141

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	180
Protein 2 (blue)	135
Protein 3 (blue)	100
Protein 4 (orange)	75
Protein 5 (blue)	63
Protein 6 (blue)	48
Protein 7 (blue)	35
Protein 8 (green)	25
Protein 9 (blue)	17
Protein 10 (blue)	11
Storage (°C)	-15 to -25
Minimum Shelf Life	05.06.2026

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.02.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.