

**CERTIFICATE OF ANALYSIS**

<b>SERVA Triple Color Protein Standard III</b>	<b>Cat.No. : 39258</b>
5 – 245 kDa, Liquid Mix	<b>Lot No. : 161186</b>

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	245
<b>Protein 2 (blue)</b>	180
<b>Protein 3 (blue)</b>	135
<b>Protein 4 (blue)</b>	100
<b>Protein 5 (orange)</b>	75
<b>Protein 6 (blue)</b>	63
<b>Protein 7 (blue)</b>	48
<b>Protein 8 (blue)</b>	35
<b>Protein 9 (green)</b>	25
<b>Protein 10 (blue)</b>	20
<b>Protein 11 (blue)</b>	17
<b>Protein 12 (blue)</b>	11
<b>Protein 13 (blue)</b>	5
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	23.12.2018

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 04.01.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.