

**CERTIFICATE OF ANALYSIS**

---

<b>SERVA Triple Color Protein Standard III</b>	<b>Cat.No. : 39258</b>
5 – 245 kDa, Liquid Mix	<b>Lot No. : 200188</b>

---

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	245
<b>Protein 2 (blue)</b>	180
<b>Protein 3 (blue)</b>	135
<b>Protein 4 (blue)</b>	100
<b>Protein 5 (orange)</b>	75
<b>Protein 6 (blue)</b>	63
<b>Protein 7 (blue)</b>	48
<b>Protein 8 (blue)</b>	35
<b>Protein 9 (green)</b>	25
<b>Protein 10 (blue)</b>	20
<b>Protein 11 (blue)</b>	17
<b>Protein 12 (blue)</b>	11
<b>Protein 13 (blue)</b>	5
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	20.02.2022

<b>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</b>
--

---

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date: 21.02.2020**

Christian Monsler

Marina Konter

This report has been computer generated and does not contain a signature.