

CERTIFICATE OF ANALYSIS

SERVA Triple Color Protein Standard III 5 – 245 kDa, Liquid Mix	Cat.No. : 39258 Lot No. : 201033
---	---

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	245
Protein 2 (blue)	180
Protein 3 (blue)	135
Protein 4 (blue)	100
Protein 5 (orange)	75
Protein 6 (blue)	63
Protein 7 (blue)	48
Protein 8 (blue)	35
Protein 9 (green)	25
Protein 10 (blue)	20
Protein 11 (blue)	17
Protein 12 (blue)	11
Protein 13 (blue)	5
Storage (°C)	-15 to -25
Minimum Shelf Life	02.10.2022

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 06.10.2020

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.