

CERTIFICATE OF ANALYSIS**SERVA Pink Color Protein Standard II**
10.5 – 175 kDa, Liquid Mix**Cat.No. : 39259**
Lot No. : 150047

Protein Components	Molecular Weight (kDa)
Protein 1 (pink)	175
Protein 2 (pink)	130
Protein 3 (blue)	95
Protein 4 (pink)	70
Protein 5 (pink)	62
Protein 6 (pink)	51
Protein 7 (blue)	42
Protein 8 (pink)	29
Protein 9 (pink)	22
Protein 10 (pink)	14
Protein 11 (blue)	10.5
Storage (°C)	-15 to -25
Minimum Shelf Life	12.01.2018

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control**Printing date: 25.10.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.