

CERTIFICATE OF ANALYSIS

| | |
|---|---|
| SERVA Pink Color Protein Standard II 10.5 – 175 kDa, Liquid Mix | Cat.No. : 39259 Lot No. : 170443 |
|---|---|

| Protein Components | Molecular Weight (kDa) |
|---------------------------|-------------------------------|
| Protein 1 (pink) | 175 |
| Protein 2 (pink) | 130 |
| Protein 3 (blue) | 95 |
| Protein 4 (pink) | 70 |
| Protein 5 (pink) | 62 |
| Protein 6 (pink) | 51 |
| Protein 7 (blue) | 42 |
| Protein 8 (pink) | 29 |
| Protein 9 (pink) | 22 |
| Protein 10 (pink) | 14 |
| Protein 11 (blue) | 10.5 |
| Storage (°C) | -15 to -25 |
| Minimum Shelf Life | 05.05.2020 |

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 08.04.2019

Christian Monsler

Marina Konter

This report has been computer generated and does not contain a signature.