

CERTIFICATE OF ANALYSIS

SERVA VisiBlot Standard I	Cat.No. : 39260
25 – 150 kDa, Liquid Mix	Lot No. : 160047

Protein Components	Molecular Weight (kDa)
Protein 1	150
Protein 2	120
Protein 3	90
Protein 4 (orange)	85
Protein 5	75
Protein 6	60
Protein 7 (orange)	45
Protein 8	40
Protein 9 (blue)	25
Protein 10	25
Storage (°C)	-15 to -25
Minimum Shelf Life	21.12.2018

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 26.09.2017

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.