

# CERTIFICATE OF ANALYSIS

**SERVA VisiBlot Standard I**  
25 – 150 kDa, Liquid Mix

**Cat.No. : 39260**  
**Lot No. : 230880**

| <b>Protein Components</b> | <b>Molecular Weight<br/>(kDa)</b> |
|---------------------------|-----------------------------------|
| <b>Protein 1</b>          | 150                               |
| <b>Protein 2</b>          | 120                               |
| <b>Protein 3</b>          | 90                                |
| <b>Protein 4 (orange)</b> | 85                                |
| <b>Protein 5</b>          | 75                                |
| <b>Protein 6</b>          | 60                                |
| <b>Protein 7 (orange)</b> | 45                                |
| <b>Protein 8</b>          | 40                                |
| <b>Protein 9 (blue)</b>   | 25                                |
| <b>Protein 10</b>         | 25                                |
| <b>Storage (°C)</b>       | -15 to -25                        |
| <b>Minimum Shelf Life</b> | 07.02.2026                        |

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 25.07.2025**

Dr. Judith Koch

Marina Konter

This report has been computer generated and does not contain a signature.