

**CERTIFICATE OF ANALYSIS**

<b>SERVA Triple Color Protein Standard I</b>	<b>Cat.No. : 39262</b>
15 – 180 kDa, Liquid Mix	<b>Lot No. : 230684</b>

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	170
<b>Protein 2 (blue)</b>	130
<b>Protein 3 (blue)</b>	100
<b>Protein 4 (red)</b>	70
<b>Protein 5 (peacock green)</b>	60
<b>Protein 6 (blue)</b>	45
<b>Protein 7 (blue)</b>	35
<b>Protein 8 (blue)</b>	25
<b>Protein 9 (green)</b>	15
<b>Storage (° C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	11.10.2025

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 21.06.2023**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.