

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

|                                         |                          |               |
|-----------------------------------------|--------------------------|---------------|
| <b>SERVA FastLoad 100 bp DNA Leiter</b> | <b>Kat.Nr./Cat.No. :</b> | <b>39316</b>  |
| SERVA FastLoad 100 bp DNA ladder        | <b>Lot/Contr.No. :</b>   | <b>200303</b> |

| <b>Parameter</b><br>parameter                   | <b>Methode</b><br>method                                                                                                                                                                                                  | <b>Spezifikation</b><br>specification | <b>Ergebnis</b><br>result |
|-------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|---------------------------|
| <b>Beschreibung</b><br>description              | Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 12 Fragmenten im Bereich 100 Bp bis 3000 Bp.<br>Ready-to-Use DNA standard for agarose gel electrophoresis containing 12 fragments from 100 – 3000 base pairs. |                                       |                           |
| <b>Anwendungstest</b><br>application test       | Elektrophorese<br>Electrophoresis                                                                                                                                                                                         | entspricht<br>corresponds             | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life |                                                                                                                                                                                                                           |                                       | 23.03.2022                |
| <b>Lagerung (°C)</b><br>storage                 |                                                                                                                                                                                                                           |                                       | -15 bis/ to -25           |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 24.03.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.  
This report has been computer generated and does not contain a signature.