

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

SERVA FastLoad 1 kb DNA Leiter	Kat.Nr./Cat.No.:	39317
SERVA FastLoad 1 kb DNA ladder	Lot/Contr.No. :	180156

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 14 Fragmenten im Bereich 250 Bp bis 25 000 Bp. Ready-to-Use DNA standard for agarose gel electrophoresis containing 14 fragments from 250 – 25000 base pairs.		
Anwendungstest application test	Elektrophorese Electrophoresis	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			05.02.2020
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.02.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.