

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>SERVA FastLoad 1 kb DNA Leiter</b>	<b>Kat.Nr./Cat.No. :</b>	<b>39317</b>
SERVA FastLoad 1 kb DNA ladder	<b>Lot/Contr.No. :</b>	<b>210697</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 14 Fragmenten im Bereich 250 Bp bis 25 000 Bp. Ready-to-Use DNA standard for agarose gel electrophoresis containing 14 fragments from 250 – 25000 base pairs.		
<b>Anwendungstest</b> application test	Elektrophorese Electrophoresis	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			31.05.2024
<b>Lagerung (°C)</b> storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 15.06.2021

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.