

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA FastLoad 1 kb DNA Leiter	Kat.Nr./Cat.No. :	39317
SERVA FastLoad 1 kb DNA ladder	Lot/Contr.No. :	221072

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 14 Fragmenten im Bereich 250 Bp bis 25 000 Bp. Ready-to-Use DNA standard for agarose gel electrophoresis containing 14 fragments from 250 – 25000 base pairs.		
Anwendungstest application test	Elektrophorese Electrophoresis	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			07.11.2024
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.11.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.